TOC

Thor Lichtenthaler
Produktspesialist for spektroskopi og elementanalyse
Matriks 2018

Rent vann er en begrenset ressurs

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
BOF – Biologisk oksygen forbruk (BOD)

1. Bacteria

Oxygen saturated bacteria solution is added to the sample

2.

Determination of the oxygen contents

3. T = const.

The closed sample is kept at 20°C. A part of the organic compounds are metabolized by bacteria under consumption of dissolved oxygen.

4.

After 5 days the oxygen concentration is determined. The BOD_5 in mg/L O_2 is calculated from the consumed oxygen.

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


KOF – Kjemisk Oksygen Forbruk (COD)

1.  Sulfuric acid, silver ions (catalyst) and potassium dichromate (oxydation agent) are added to the sample.

2.  The solution is refluxed at 148°C (boiling point) for 2 hours to oxidized dissolved organic compounds.

$\text{organic substance} + \text{Cr}_2\text{O}_7^{2-} \rightarrow \text{CO}_2, \text{H}_2\text{O}, \text{Cr}^{3+}$

3.  After cooling unconsumed Dichromate is determined by back titration with Fe^{2+} . The chemical oxygen demand is calculated from the difference with the consumed oxydation agent.

DIN 38 409 H41

Disadvantage:

Reagents:

concentrated sulfuric acid H_2SO_4 – corrosive!

potassium dichromate $\text{K}_2\text{Cr}_2\text{O}_7$ - toxic!

mercury substance (masking of chloride) – toxic!

Analysis time: > 2 hours, automation very difficult

Påstand (fra TOC leverandører)

TOC means *Total Organic Carbon*. The parameter TOC gains more and more importance compared to the *Chemical Oxygen Demand* (COD) since COD can only be automated with high expenditures. Furthermore, large amounts of aggressive, respectively toxic chemicals are necessary. While measuring the COD all other oxidizable compounds will also be determined in one sample – also those which are not of organic nature.

TOC som alternativ

U.S. Safe Drinking Water Act in 2001, TOC analysis emerged as a quick and accurate alternative to the classical but more lengthy biological oxygen demand (BOD) and chemical oxygen demand (COD) tests traditionally reserved for assessing the pollution potential of drinking and wastewaters.

Safe Drinking Water Act

Long title	An Act to amend the Public Health Service Act to assure that the public is provided with safe drinking water, and for other purposes
Nicknames	SDWA
Enacted by	the 93rd United States Congress
Effective	December 16, 1974

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Total inorganic carbon (TIC)

Total organic carbon (TOC)

POC purgeable organic carbon (also called **VOC** - volatile organic carbon)

NPOC non purgeable organic carbon (also called **NVOC** – non volatile organic carbon)

DOC dissolved organic carbon - organic carbon that is completely in solution. This carbon can be measured like NPOC after filtration of the sample through a 45 µm mesh. The remainder is classified as **NDOC** (non dissolved organic carbon).


BDOC biodegradable organic carbon - organic carbon that may be metabolized by bacteria or other lifeforms

AOC assimilable organic carbon - organic carbon that may be incorporated by bacteria or other lifeforms

Total Organic Carbon

TC [total carbon] = **TIC** [total inorganic carbon] + **TOC** [total organic carbon]

TOC = **POC** [purgable organic carbon] + **NPOC** [non-purgable organic carbon]

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$$\text{POC} \cong 0 \Rightarrow \text{TOC} = \text{NPOC}$$

$$\text{or } \text{TC} = \text{NPOC} + \text{TIC}$$

$$\text{Differential mode: } \text{TOC} = \text{TC} - \text{TIC}$$

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Method	Advantages	Disadvantages
Differential:	Includes volatile organic compounds	Less suitable for high TIC concentrations (TOC << TIC) long measuring times
Direct: NPOC	Suitable for high TIC concentrations (TIC >> TOC) fast results	Does not include volatile organic compounds if present in the sample

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Drinking water:	0.2 – 5 ppm TOC
Surface water:	0.5 – 50 ppm TOC
Waste water:	5 – 10,000 ppm TOC
Water in pharmacy:	0.01 – 1 ppm TOC
Water in power plants:	0.01 – 10 ppm TOC
Water in semiconductor industry:	0.001 – 0.5 ppm TOC
Waste and soil (solid TOC):	0.01 to 50% TOC

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



Ikke for faste stoffer

IR deteksjon

**New wide range IR
detector:
Low ppb to % levels**

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

Matriks uavhengig analyse

Svært stabil kalibrering over tid - dekker ett stort dynamisk måleområde

Longterm stability and reliability of the calibration

Stability of the calibration is in the range of several months up to years. The straightforward and robust basic concept make the ELEMENTAR analysers real working horses in the laboratory routine.

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
PARAMETER	DESCRIPTION
TOC ₄₀₀	total organic carbon which is released up to 400°C
ROC *	residual oxidizable carbon, determined between the signal minimum at 400°C and at 600°C in the case of dry combustion in a current of oxygen
TIC ₉₀₀	total inorganic carbon which is released up to 900°C
TOC	total organic carbon, sum of TOC ₄₀₀ and ROC
C	total carbon, sum of TOC ₄₀₀ , ROC and TIC ₉₀₀

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soil TOC cube – separation of carbon fractions

3-step temperature ramp

Soil samples (Cambisol B Horizon, 150 mg each)

New parameters for the determination of carbon

TOC₄₀₀: **T**otal **O**rganic **C**arbon (released between 150 and 400° C in the presence of O₂)

ROC: **R**esidual **O**xidisable **C**arbon (elemental carbon, e.g., carbon black, soot, released between 400° C and 600° C in the presence of O₂)

TIC₉₀₀: **T**otal **I**norganic **C**arbon (released between 600 and 900° C in the presence of O₂)

Parameter	Start Temp. [° C]	Heating Rate [° C/min.]	Target Temp. [° C]	Hold time [s]
TOC ₄₀₀	150	70	400	120
ROC	400	70	600	100
TIC ₉₀₀	600	70	900	100

Significant cost savings for waste disposal on the horizon...

...when elemental carbon is substracted from TOC

SUBSTANCE	TOC ₄₀₀ [%]	ROC [%]	TIC ₉₀₀ [%]
FLUVISOL	3.63	0.25	0.041
ARENOSOL	1.65	0.087	0.008
SOIL WITH SLUDGE	2.06	0.12	0.033
EXCAVATED SOIL	0.23	0.15	0.42
COAL MINE TAILINGS	17.21	18.93	11.55
FOUNDRY SAND	1.21	1.93	1.27
WASTE INCINERATOR ASH	0.51	0.52	0.19
SLAG	0.024	0.098	0.16
CONTROL MIXTURE ACCORDING TO DIN 19539	2.14	1.73	2.08


- Easy to decompose: Calcite (10% HCl, RT)

- Difficult to decompose: Magnesite, Dolomite, Siderite (10% HCl, 50° C). During heating also organic material can be decomposed to form volatile compounds.

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
EXCELLENCE IN ELEMENTS

www.elementar.de

110 years tradition... Heraeus

Wilhelm Carl Heraeus

Elementar´s TOC product range

vario TOC select / cube

- high temperature combustion technique
- liquid samples
- solid samples
- TC, TOC, TIC, POC
- TN with additional detectors, EC, CLD, IR

acquray

- UV / persulfate oxidation
- liquid samples
- solid samples with additional solids module
- TC, TOC, TOC, TIC
- TN and TP with additional modules

Soli TOC cube

- high temperature combustion technique
- temperature ramp method
- solid samples
- TC, TOC, TIC, ROC
- TN with additional EC detector

2016 – Nye lokaler og Fabrikk i Langenselboldt nær Hanau - Tyskland

Elementar moves into new headquarters
in Langenselboldt near Frankfurt, Germany

Matriks samarbeid siden 2009

Good to know,
that we are a strong family.

As a family enterprise we stand for
highest reliability and personal commitment.